



M1-540-15FM

Single-Stage Preset Flowmeter Regulator with CGA-540 Nut and Nipple Inlet



NMR-870-15FM

Single-Stage Flowmeter Regulator with CGA-870 Yoke

OXYGEN MR CONDITIONAL SLIM-LINE FLOWMETER

MODEL	FLOW RANGE	REGULATOR PRESSURE	INLET
NMR-870-15FM	1/2 - 15 LPM	50 PSI preset	CGA-870 Yoke ¹

FEATURES:

- Chrome-plated brass body with all brass high-pressure chamber
- Maximum rated inlet pressure 3000 PSI
- Durable neoprene diaphragm
- Internal reseating relief valve protects against over-pressurization
- Sintered filter for additional safety and to extend regulator life
- 2" diameter gauge
- 15 LPM Oxygen models provide 75 LPM flood minimum
- 7 LPM models restricted to 10 LPM maximum
- Polycarbonate inner and outer tubes provide greater accuracy and durability
- Includes plastic swivel hose barb (not shown)
- 6 year limited warranty

SINGLE STAGE MODELS

GAS/MODEL	FLOW RANGE	INLET CONNECTION
OXYGEN		
M1-540-8FM	1/2 - 7 LPM	CGA-540 Nut and Nipple ¹
M1-540-8FMH	1/2 - 7 LPM	CGA-540 Handtight Nut and Nipple ¹
M1-870-8FM	1/2 - 7 LPM	CGA-870 Yoke ¹
M1-540-15FM	1/2 - 15 LPM	CGA-540 Nut and Nipple ¹
M1-540-15FMH	1/2 - 15 LPM	CGA-540 Handtight Nut and Nipple ¹
M1-870-15FM	1/2 - 15 LPM	CGA-870 Yoke ¹
M1-540-15FMPTO	1-15 LPM	CGA-540 Nut and Nipple w/Power Take Off ¹
MEDICAL AIR		
M1-346-15FM	1 - 15 LPM	CGA-346 Nut and Nipple ¹
M1-346-15FMH	1 - 15 LPM	CGA-346 Handtight Nut and Nipple ¹
M1-950-15FM	1 - 15 LPM	CGA-950 Yoke ¹
CARBON DIOXIDE		
M1-320-12FM	1/2 - 12 LPM	CGA-320 Nut and Nipple ¹
M1-320-12FMH	1/2 - 12 LPM	CGA-320 Handtight Nut and Nipple ¹
M1-940-12FM	1/2 - 12 LPM	CGA-940 Yoke ¹

See full CGA chart on inside back cover

FLOWMETER FEATURES:

- Polycarbonate inner and outer tubes
- Calibrated at 50 PSI inlet
- DISS 1240 outlet

REGULATOR FEATURES:

- All-Brass body and high-pressure chamber
- Internal reseating relief valve protects against over-pressurization
- 2" diameter gauge
- Sintered filter for additional safety and extended regulator life
- Durable neoprene diaphragm
- 3000 PSI maximum inlet pressure
- DISS 1240 outlet
- 50 PSI preset outlet

¹ **⚠ WARNING:** This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.